

Return & Repair Form:

Customer Information:

Company Name: Delanco Twp B: _____ S: _____

Date Received: 12-14-09 Date Given to Service: _____

Carrier: FedEx UPS - DHL - USPS Method: GRD NDA 2DY - 3DY - Other

Product Information:

Product: 7410-7110-6510-6810-7510
9510- Drug Tester 5000

Priority
Serial #: AR WE-0015

Description: A - B - Plus - Screener - Demo

Printer Ser #: AR -

Whole Instrument

Top ½

Sim Ser #: _____

Other _____

Probe Ser # DD P -

ACCESSORIES:

☐ 110 V A/C Adapter ☐ Regulator ☐ Mag Card Rdr# _____

☐ Printer Paper ☐ Printer Ribbon ☐ Casio # _____

☐ Mouthpieces ☐ Carrying Case ☐ Dry Gas

Other (Please Specify) _____ Warranty Exp. Date _____

Repair Information:

Test #

Reason for Return:

Part Number	Description	Qty	Total Cost
<u>12029</u>	<u>Power Supply</u>	<u>1</u>	<u>N/A</u>
<u>6808455</u>	<u>FC</u>	<u>1</u>	<u>N/A</u>
<u>6808486</u>	<u>Plate</u>	<u>1</u>	<u>N/A</u>
<u>MP Cal 71</u>	<u>Cal</u>	<u>1</u>	<u>N/A</u>
<u>MP Labor</u>	<u>Labor</u>	<u>.5</u>	<u>N/A</u>

Repair Notes:

Unit arrived with bad Power supply, Install new power supply.

FC ARWL-0217 out of date, Install FC ARHJ-0169

Cal w/QC; Ops Check

Service Technician: JS

Date: 12/15/09